


UlceRx Benefits:

- **Heal faster:** *UlceRx is scientifically proven to heal venous ulcers of 4cm or less in half the time than with traditional bandages (ulcers >4cm healed in the same time as with bandages).*
- **Reduce pain:** *patients using the UlceRx reported no pain at night or in the morning.*
- **Easy wear:** *no complicated bandages, just a two-step process to put on.*

ULCER_X

Therapy Solution



SIGVARIS, Inc.
1119 Highway 74 South
Peachtree City, GA 30269
USA
1.800.322.7744

SIGVARIS Corp.
4535 Dobrin
Ville St-Laurent, Quebec
H4R 2L8
Canada
1.800.363.4999

Study cited: The treatment of venous leg ulcers with a specially designed compression stocking kit. F. Mariani, et al. *Phlebologie* 2008; 37: 191-197. September 2008.

for more information,
please visit
www.sigvarisusa.com
www.sigvaris.ca

SIGVARIS and LIFE FOR LEGS are registered trademarks of Ganzoni & Cie AG,
CH-9014 St.Gallen/Switzerland, in many countries worldwide

LIT-UXBROC REV. 06.09

SIGVARIS

LIFE FOR LEGS

Venous ulcers
heal faster and
with less pain
with **ULCER_X**
Therapy Solution

Proven superior
treatment choice
for the successful
healing of
venous ulcers

The SIGVARIS UlceRx: *designed with patients in mind*

This dual-stocking system consists of a low-compression underliner and a 30-40 mmHg overstocking. The underliner keeps your dressing in place at night, helping you rest comfortably. Its smooth exterior also aids in the donning of the overstocking. During your daily activities, the overstocking provides the necessary compression needed to heal your venous ulcer.

Scientifically Proven Results

- **Complete wound closure achievement:**
 - » 96.2% with UlceRx Kit as compared to 70% with traditional bandages
- **Ulcers with diameter of up to approximately 4cm healed twice as fast with the UlceRx Kit compared to bandages**
 - » Larger ulcers (>4cm) healed as rapidly with the UlceRx Kit as with bandages
- **Reported pain at night/morning was absent with the UlceRx Kit group but reported by 40%/20%, respectively, in the bandage group**

The UlceRx Kit

The UlceRx Kit is made of two easy-to-use components: a low-compression underliner and a 30-40 mmHg overstocking. The underliner is made with an inner layer of cotton and a smooth exterior that allows easy donning of the overstocking. The 30-40 mmHg overstocking is recommended for removal at night, but you should wear the low-compression understocking during sleep or while recumbent.



Underliner

+



503 Overstocking

= **ULCER_x**

The SIGVARIS UlceRx makes it easy for you to successfully complete your venous leg ulcer treatment.



"Discomfort during the night and in the morning was absent with the UlceRx Kit."

The treatment of venous leg ulcers with a specially designed compression stocking kit. F. Mariani, et al. Phlebologie 2008; 37: 191-197. September 2008.